







Same open length for FMS & FWS series	No. of panels	Part no.	Length (When open)		Height (w/ 3" casters)		Folded size		Weight	
			cm	inch	cm	inch	Width cm	Depth cm	FMS	FWS
									kg	
	4	FMS14504 FWS14504	120	47.2"	145	57"	30	7.2	23.6	11.2
		FMS16504 FWS16504			165	65"			24.7	12.5
		FMS18504 FWS18504			185	73"			26.1	13.9
	6	FMS14506 FWS14506	180	70.9"	145	57"	30	10.8	27.2	15.5
		FMS16506 FWS16506			165	65"			29.6	17.5
		FMS18506 FWS18506			185	73"			31.5	19.3
	8	FMS14508 FWS14508	240	94.5"	145	57"	30	14.4	32.0	19.9
		FMS16508 FWS16508			165	65"			34.6	22.4
		FMS18508 FWS18508			185	73"			36.9	24.8
	10	FMS14510 FWS14510	300	118.1"	145	57"	30	18	35.8	23.6
		FMS16510 FWS16510			165	65"			38.9	26.7
		FMS18510 FWS18510			185	73"			N.A	29.6
	12	FMS14512 FWS14512	360	141.7"	145	57"	30	21.6	40.1	28.1

screen color options:




Machan International Co., Ltd.

No.354, Dama Rd., Waipu Dist.,
Taichung City 43856, Taiwan

Tel: +886-4-2683-4666
Fax: +886-4-2683-4672

sales@bailida-medical.com
www.bailida-medical.com



2019 copyright by BAILIDA

FMS

FWS

DUAL-USE FOLDING
SCREEN

LIGHT WEIGHT
CHANGEABLE SCREEN LENGTH
OPTIMUM HYGIENE





BAILIDA



Optimum Hygiene


Traditional hospital privacy curtain is hard to clean on a frequent basis. The study shows that it can be contaminated with deadly drug-resistant bugs within **2 weeks**, including MRSA. According to the US Centers for Disease Control, nearly **20,000** MRSA-related deaths occurred in the US in 2017, most of them from hospital-acquired infections.

Bailida's folding screen is the best alternative to the traditional privacy curtain. With its **anti-bacterial** surface and quick-to-clean nature, it can efficiently minimize the risk of posing threat to patient safety.

With swivel castors attached, it provides excellent maneuverability which enables health care workers to move the folding screens around for different

applications. The elegant design is ideal for bringing a warm feeling to the hospital.

With its sturdy aluminum frame and patented structure, Bailida's folding screen is not only one of the lightest screens in the market, it also ensures the optimum hygiene. Furthermore, the panels are also impact resistant, scratch resistant and flame retardant.



• Traditional hospital privacy curtain

VS.




• Bailida's Folding Screen

Combination of lightness and panel changeable feature

The patented structure ensures:

1. Panels can be added or reduced meanwhile maintains the lightest weight.
2. No gap or light goes through panels, patient privacy is 100% guaranteed.
3. Smooth folding: 200,000 times smooth folding is guaranteed.

EX:



Quick Storage

The folding screen can be tied with an adjustable elastic band which can be used for different numbers of panels. The unique patented design ensures it does not rebound when the number of panels changed.




Dual Use : mobile and wall mounted

The screen can be detached by pushing down the click. This enables the screen to be easily transformed into its wall-mounted function.




Standing post is available for connecting with wall-mounted screen to be used as a fixed screen.



Stable Movement

The built-in shock-proof function provides additional stability over uneven or rough surfaces, it prevents screens from falling down during movement.



Quick To Clean

The castors can be removed for cleaning purpose.

